

AllSet+™ Gold SSP
HLA-ABDRDQ WORKSHEET
Code no.: 54360D

Name _____

Sample I.D. _____

Test Date _____

Serological Equivalent	Allele Specificity	24	25	26	27	28	29	30
DQ2, -	DQB1*02:01:01-11			26				
DQ7, DQ8, -	DQB1*03:01:01:01-01-01:08/02:05/04/09-10/13-14/16/19/21-24/27-29/35-36/42/44/46-47				27			
DQ8, -	DQB1*03:02:01-02:04/05:01-05:04/07-08/11/18/32/37/45					28		
DQ9, DQ3, -	DQB1*03:03:02:01-03:04/06/12/15/17/20/25-26/30-31/33-34/38-41/43, DQB1*04:03:02						29	
DQ4, -	DQB1*04:01:01-02:01/03:01/04-09						29	30
DQ5, -	DQB1*05:01:01-01-11/13-18	24						
DQ6, -, DQ1, null	DQB1*06:01:01-50/52		25					
-	DQB1*06:51		25				29	
Product sizes			230	205/250	195	140/205	70	120/125/140/145/150
							205	
	Lane Number	24	25	26	27	28	29	30
	Positive Lanes							
	1st match:							
	2nd match:							
	Failed Controls							
	False Positive							
	False Negative							

AllSet™ Gold SSP
HLA-ABDRDQ WORKSHEET
Code no.: 54360D

Name _____

Sample I.D. _____

Test Date _____

Serological Equivalent	Allele Specificity	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53
A11, -, null	A*11:01:01-01:20/01:22-01:25/01:27-01:43/01:45-03/05/07-16/20-25/29-34/36-37/41-42/44-49/51-56/58-89/91-93/95-97/99N-105/107-111/114-117/119-120/122-129/131-132/135-138				34														48	49				
A11, -	A*11:01:21/01:44/04/19/27/35/38-39/90/98/106/121				34															49				
A11, -	A*11:01:26/06				34														48					
-	A*11:17/40				34				39															
-	A*11:26/118, A*66:18																		48	49				
-	A*11:43				34						42								48	49				
-, A66(34), A10, A26	A*11:50Q/94, A*66:02-03/05/10/15-16					36														49				
-	A*11:57				34														48	49	50			
-	A*11:112				34	36													48	49				
-	A*11:113				34						41								48	49				
-	A*11:133				34														48					
-	A*11:139								37										48	49		51		
A23, -, null	A*23:01:01-03:02/05-12/15-33/35-44/46-50/52-56					35																		
A23	A*23:04					35										45								
-	A*23:13					35																51		
-	A*23:14:01-14:02, A*24:71					35	37																	
-	A*23:45					35												47						
A24, -, A2403, null	A*24:02:01:01-02:25/02:27-02:36/02:38-02:64/04-06/08/10:01-17/20-21:02/23/25-32/34-43/45N-50/52-61/63-70/72-81/83N-91/93/95-107/109-111/113-124/126-130/132N-137/139-166/168-187/189-203/205-206/208-212/214-227/229-230						37																	
A24	A*24:02:26						37																	
A24	A*24:02:37						37																	
A2403, A9, -	A*24:03:01-03:02/22/33/125/167/207						37								45									
A24, -	A*24:07/108/112/131						37																51	
null	A*24:09N						37										46							
-	A*24:18/204						37									45	47							
A9	A*24:19						37													49		51		
-	A*24:24																					51		
-	A*24:44						37														49			
-	A*24:51/92				34		37																	
-	A*24:62/213						37											47						
-	A*24:82						37					42	44											
-	A*24:94/138															45						50		
A25, -, null	A*25:01:01-01:07/03-04/07-22							38																
A25	A*25:02							38												48				
-	A*25:05							38														50		
-	A*25:06							38								45								
A26, A10, null, -	A*26:01:01:01-03:01:15/01:17-01:22/01:24-02/04/07:01-08/10-12/14-18/20/22-29/31-32/34-53/55-77/79-84								39										47					
A26	A*26:01:02/05								39															
A26	A*26:01:16								39										47					
A26	A*26:01:23					36		39											47					
A26, -	A*26:03:01-03:02/06/21/30/78					36																		
A26	A*26:09								39							45			47					
-	A*26:13								39										47	48				
Product sizes					155/195/205/25/565/575		225/307/35																	
Lane Number		31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53
Positive Lanes																								
1st match:																								
2nd match:																								
Failed Controls																								
False Positive																								
False Negative																								

**AllSet+™ Gold SSP
HLA-ABDRDQ WORKSHEET
Code no.: 54360D**

Name _____

Sample I.D. _____

Test Date _____

Serological Equivalent	Allele Specificity	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53																						
A68, -, A28, null	A*68:01:01:01-01:10/01:12-04/06-14/16-19/21:01-28/30-35/37-47/51-60/62-70/72-82/85-99																				49	50																								
-	A*68:29														44						49	50																								
-	A*68:61																			48	49	50																								
-	A*68:71			33																	49	50																								
-	A*68:84									39												50																								
A69	A*69:01																				49	51																								
A74, -, null	A*74:01-05/07-09/11-12N/14N-17																				49			52																						
A74, C6	A*74:06, C*06:02:15																							52																						
-	A*74:10									40									47					52																						
-	A*74:13			33																	49			52																						
A80, -	A*80:01:01:01-02																							53																						
Product sizes				155/195/205/525/565/575 640		229/230/735		170/215/230		530/535		490/520/550/580		170/210/510/550		385/440		75/235		115/150		420/530/565/675		200/205/250		440		170/180/205/210		170/570		325/640		385/400/420/425/445/460/525/560/565		550		180		155/180		350/380/385		155/180		495
Lane Number		31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53																						
Positive Lanes																																														
1st match:																																														
2nd match:																																														
Failed Controls																																														
False Positive																																														
False Negative																																														

